

LABORATORY REPORT

RIQAS®

PATIENT NAME : MR K G PREMASIRI

AGE : 81 Y 0 M 0 D BILL NO : 255245 (Cash [Full])

GENDER : Male REPORT DATE & TIME : Thu 15 Oct, 2026 09:37 PM

SPECIMEN : Serum SAMPLE DATE & TIME : Thu 15 Oct, 2026 08:16 PM

REFERRED BY : DR CHINTHAKA DE SILVA (CC) |

CENTER : COOP-HOSPITAL MATARA

TROPONIN I TITRE (QUANTITATIVE)

TEST	RESULT	UNIT	REF. RANGE
TROPONINE I - TITRE	0.030	ng/mL	

Reference Ranges

0.0 - 0.06ng/ml - Negative

0.06- 0.1 ng/ml - High risk; advisable to repeat the test within 1hr with fresh specimen

Above 0.1 ng/ml - Positive

Comments :

detects lower levels of Troponin I, than conventional Troponin I assays, hence earlier detection of onset.

Troponin I above 0.06 (99th percentile URL) Indicates myocardial necrosis. Troponin I <0.06 does not exclude a clinical diagnosis of acute coronary syndrome (myocardial damage). Suggest repeat quantitative hS Troponin I after 1 to 2 hours to exclude an ongoing myocardial necrosis in highly suspicious cases.

Third Universal Definition of Myocardial infarction, circulation. 2012; 126

Analytical sensitivity : 0.06 ng/mL

Enzyme Immunoassay method specific reference ranges given.

Test performed using TOSOH AIA - 360 Fully Automated Immunoassay Analyzer

Dr. Thushari Vithanage
MBBS, PG Dip. Path & MD (Chemical Pathology)
Consultant Chemical Pathologist
DGH Matara

Dr. Dileepa Amarasekiriya
MBBS, D. Path, MD (Haematology)
Consultant Haematologist
DGH Matara

Dr. N J Pittigalage
MBBS, Dip. Microbiology, MD. Microbiology
Consultant Microbiologist
DGH Matara



Medical Laboratory Technologist

No 223, Kumarathunga Mawatha, Matara
printed by: 69 (Thu 2026-10-15 09:42pm)

Tel : +94(0)412232893

Manage by MarCoop health care products private Ltd.

Email: cooplabmatara@gmail.com

Web: www.cooperativehospitals.lk